

ACCOMMODATIVE ORTHOTICS

TYPE	MATERIALS	PROPERTIES
1	EVA 45 EVA 65	<ul style="list-style-type: none"> • THE MOST RIGID • A GOOD COMPROMISE TO THE FUNCTIONAL ORTHOTIC
2	EVA 35 EVA 65	<ul style="list-style-type: none"> • FOR MINOR DEFORMITY TO MODERATE (MINOR RHUMATOID ARTHRITIS, NON PATHOLOGICAL DIABETIC FOOT, LITTLE NEOPLASM...)
3	EVA 35 EVA 35	<ul style="list-style-type: none"> • FOR MINOR DEFORMITY TO MODERATE (RHUMATOID ARTHRITIS WITH DIGITAL DEVIATIONS, NON PATHOLOGICAL DIABETIC FOOT, LITTLE NEOPLASM...)
4	P-CELL 15 EVA 45	<ul style="list-style-type: none"> • FOR MODERATE DEFORMITY (RHUMATOID ARTHRITIS WITH digital DEVIATIONS AND ATROPHY OF THE FAT PAD, DIABETIC FOOT WITH ULCER RISK...)
5	P-CELL 15 EVA 35	<ul style="list-style-type: none"> • FOR MODERATE DEFORMITY TO SEVERE (RHUMATOID ARTHRITIS WITH DEVIATIONS AND ATROPHY OF THE FAT PAD, DIABETIC FOOT WITH ULCER, CHARCOT FOOT...)
6	PLASTAZOTE P-CELL 15	<ul style="list-style-type: none"> • FOR SEVERE DEFORMITY OR HYPERSENSITIVE SKIN OR BIGGER NEOPLASM.

PodoKop inc.

17805 Lapointe

Mirabel, Qc, J7J 0W7

☎ (450) 979-1772 ✉ info@podokop.com

